# ***Mode d'emploi***

---

*MNPG93 Rev. 01 Edition 05/05/15*

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*LymphoPress 4*

**I.A.C.E.R. Srl**  
**Via S. Pertini 24/A – 30030**  
**Martellago (VE) ITALY**  
**Tel. 041.5401356 • Fax 041.5402684**  
**e-mail: [iacer@iacer.it](mailto:iacer@iacer.it) <http://www.iacer.it>**

---

## Sommaire

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## EMC tables

<b>Guidance and MANUFACTURER'S declaration – ELECTROMAGNETIC EMISSIONS – FOR ALL THE EQUIPMENT AND THE SYSTEMS</b>		
LYMPHOPRESS 4 is intended for use in the electromagnetic environment specified below. The customer or the user of the LYMPHOPRESS 4 should assure that it is used in such an environment.		
<b>Emission test</b>	<b>Compliance</b>	<b>Electromagnetic environment – guidance</b>
RF emissions Cispr 11	Group 1	LYMPHOPRESS 4 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment
RF emissions Cispr 11	Class B	LYMPHOPRESS 4 is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A Complies	The LYMPHOPRESS 4 is suitable for use in all establishments, other than domestic establishments and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

<b>Guidance and MANUFACTURER'S declaration – ELECTROMAGNETIC IMMUNITY – FOR ALL THE EQUIPMENT AND THE SYSTEMS</b>			
LYMPHOPRESS 4 is intended for use in the electromagnetic environment specified below. The customer or the user of LYMPHOPRESS 4 should assure that is used in such environment.			
<b>Immunity test</b>	<b>Test level EN 60601-1-2</b>	<b>Compliance level</b>	<b>Electromagnetic environment – guidance</b>
Electrostatic discharge (ESD) EN 61000-4-2	± 6kV contact ± 8kV air	± 6kV contact ± 8kV air	Floors should be wood, concrete or ceramic tile. If floor are covered with synthetic material, the relative humidity should be at least 30%

Electrical transient/burst IEC 61000-4-4	fast $\pm 2\text{kV}$ for power supply lines	$\pm 2\text{kV}$ per power supply lines	Mains power quality should be at that of a typical commercial or hospital environment.
Impulses EN 61000-4-5	$\pm 1\text{kV}$ differential mode	$\pm 1\text{kV}$ differential mode	Mains power quality should be at that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	$< 5\% U_T$ (>95% dips of $U_T$ ) per 0,5 cycles  $40\% U_T$ (60% dips of $U_T$ ) per 5 cycles  $70\% U_T$ (30% dips of $U_T$ ) per 25 cycles  $< 5\% U_T$ (>95% dips of $U_T$ ) per 5 seconds	$< 5\% U_T$ (>95% dips of $U_T$ ) per 0,5 cycles  $40\% U_T$ (60% dips of $U_T$ ) per 5 cycles  $70\% U_T$ (30% dips of $U_T$ ) per 25 cycles  $< 5\% U_T$ (>95% dips of $U_T$ ) per 5 seconds	Mains power quality should be at that of a typical commercial or hospital environment.  If the user of the LYMPHOPRESS 4 requires continued operation during power mains interruptions, it is recommended that LYMPHOPRESS 4 be powered from an uninterruptible power supply or a battery.
Mains power electromagnetic field EN 61000-4-8	3 A/m	3 A/m	Mains power quality should be at that of a typical commercial or hospital environment.

**Guidance and MANUFACTURER'S declaration – ELECTROMAGNETIC IMMUNITY –  
FOR ALL THE EQUIPMENT AND THE SYSTEMS THAT ARE NOT LIFE-SUPPORTING**

LYMPHOPRESS 4 is intended for use in the electromagnetic environment specified below. The customer or the user of LYMPHOPRESS 4 should assure that is used in such environment

Immunity test	Test level EN 60601-1-2	Compliance level	Electromagnetic environment – guidance
Conducted RF EN 61000-4-6	3 Veff from 150kHz to 80MHz	3 Veff from 150kHz to 80MHz	<b>Portable and mobile RF communications equipment should be used no closer to any part of LYMPHOPRESS 4, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance:</b> $d = 1,2 \cdot \sqrt{P}$ 150kHz to 80MHz $d = 1,2 \cdot \sqrt{P}$ 80 MHz to 800 MHz $d = 2,3 \cdot \sqrt{P}$ 800 MHz to 2,5 GHz where (P) is the maximum output power rating of the transmitter in Watts (W) according to the transmitter manufacturer and (d) is the recommended separation distance in metres (m).
RF Radiata EN 61000-4-3	3 Veff from 80MHz to 2,5GHz	3 Veff from 80MHz to 2,5GHz	

Field strengths from fixed RF transmitters, are determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol:

**Recommended separation distances between portable and mobile communications equipment and the EQUIPMENT**

LYMPHOPRESS 4 is intended for the use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of LYMPHOPRESS 4 can help prevent electromagnetic interferences by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and LYMPHOPRESS 4 as recommended below, according to the maximum output power of the communication equipment.

Rated maximum power of the transmitter (W)	Separation distance according to the frequency of the transmitter (m)		
	150kHz to 80MHz $d = 1,2 \cdot \sqrt{P}$	80MHz to 800MHz $d = 1,2 \cdot \sqrt{P}$	800MHz to 2GHz $d = 2,3 \cdot \sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3

10	3,8	3,8	7,3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in Watts (W) according to the transmitter manufacturer.

Note:

- (1) At 80 MHz and 800 MHz the separation distance for the higher frequency range applies.
- (2) These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

**I.A.C.E.R S.r.l.**

*Sede operativa:*

30030 Martellago (VE) - Via. S. Pertini 24/A  
Tel +39 041 5401356 - Fax +39 041 5402684

-----  
*Sede legale:*

S. Marco 2757 - 30124 Venezia  
Cod. Fisc./P.IVA IT 00185480274  
R.E.A. VEN. 120250 - M. VE001767  
Cap.Soc. € 110.000,00 i.v.  
www.iacer.it - iacer@iacer.it

---