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
### Instructions for use

Please read before first use.  
Keep as a reference.

### Instrument identification

The name of the model and the serial number are printed on the dial of the sphygmomanometer.

Model \_\_\_\_\_


Serial number \_\_\_\_\_

### Spengler SAS

The LIAN® sphygmomanometer is a professional non-invasive medical device for measuring arterial blood pressure in conjunction with a stethoscope (auscultatory method) : the maximum systolic value or SYS (as a result of the contraction of the left ventricle) and the minimum diastolic value or DIA (as a result of relaxation of the ventricles).

### Preparations for measurement (cuff and stethoscope)

Select the correct size cuff for the patient's arm circumference (standard cuff calibration for the bladder: cuff size, index limits, artery position (middle of the bladder)).

## **Measuring:**

- wait for 5 minutes for patient's rest,
- keep the cuff at the heart level during measuring,
- keep the air release valve in close position,
- inflate the cuff to detect the estimated Systolic pressure until the radial pulse disappears (Korotkoff's sounds),
- reinflate the cuff up to 30 mmHg above estimated Systolic pressure,
- release, by opening the air release valve, at a rate of 2 mmHg per heart beat and simultaneously read the blood pressure on the manometer,
- read the SBP for phase 1 of the Korotkoff's sounds (distinct level at which the sound first appears when the air is let out),
- read the DBP for phase V of Korotkoff's sounds (where the sounds disappear again),
- measure the pulse,
- release the air fully until the complete deflation of the cuff.

**Warning: stopping bloodstream involves some risks.**

## **Hypertension**

When SBP is  $\geq 140$  mmHg and/or DBP is  $\geq 90$  mmHg.

## **Regular inspection and maintenance**

### **Metrological control by a registered association**

Check accuracy and air-leak proof:

- at least once every two years,
- after each maintenance or repair service.

### **Replacement of cuff, bladder and hand bulb:**

- check air-leak proof of bladders, tubes and hand bulbs, once a year at least,
- check the cuff appearance yearly.

### **Cleaning and disinfection:**

- manometers: wipe clean with a solution of appropriate disinfectant ; rinse with a wet non-woven tissue. Do not immerse.
- cotton or nylon cuffs: wash in machine, water not exceeding 30-35 °C. Important: the cuffs should be fastened during cleaning. Do not dry in machine.
- bladders and tubes: obturate the extremity of the tube with a clamp before washing in mild soapy water or immerse in a disinfectant solution. Rinse thoroughly and dry flat.

### **Caution**

Take care of the manometer, especially the glass (fragile).

Do not expose to heat, sun or humidity.

## Lian® sphygmomanometers

### Model

|         |  |
|---------|--|
| Classic | The twin tube "tradition" model  |
| Couleur | Robust and elegant single tube model with lateral spoon                |
| NM      | Classic single tube model: right and left hand use                     |
| Métal   | "New Look" large dial single tube model with lateral spoon             |
| Nano    | "Stylish and performing" ABS double case with spoon, single tube model |

### Bladder size range (for calibrated velcro cuff)

|                   |            |
|-------------------|------------|
| Thigh/large adult | 16 x 32 cm |
| Adult             | 12 x 27 cm |
| Child/small adult | 8 x 23 cm  |
| Infant            | 6 x 18 cm  |
| Newborn           | 4 x 12 cm  |

According to the EN CE 1060-1, the optimum width of the cuff bladder should be optimally 0.40 times the circumference of the limb at the midpoint of the intended range of the cuff. The length of the cuff bladder should be approximately 0.80 times the circumference of the limb at the midpoint of the intended range of the cuff.

### Packaging

In bag under carton sleeve. LIAN Métal under blister with bag.

### Guarantee

Spengler guarantees the products for 2 years from date of purchase (fair wear, rubber and PVC parts, cuffs are excluded). The guarantee does not apply to damages due to disassembly, shocks or misuse.

The guarantee card should be duly completed and returned to our After-Sales Service. Replacement parts are available at your Spengler distributor.

### Product quality assurance system

#### STANDARDS :

The medical devices comply with the specifications of the European norms (non-invasive sphygmomanometers EN 1060-1 and 1060-2).

#### QUALITY SYSTEM:

ISO 9001:2008

ISO 13485:2004

## CE MARK:

Aneroid sphygmomanometers certificate delivered by G-MED, registered association n° 0459. Instruments in compliance with the requirements of the Directive 93/42 EEC and its modification 2007/47/EEC.

## NF-MEDICAL MARK:

Those devices comply with the stipulations of the NF 510 certification guidelines. The main certified features are accuracy, resistance, mechanical aspects, resilience, reliability and security.

## DETAILS OF THE HOLDER:

- Identification number: 01/01
- batch number: see device's dial
- address and telephone number

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Manufacturer: Spengler S.A.S.  
28, avenue Charles-de-Gaulle  
36100 ISSOUDUN  
France

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
## Technical specifications

|                       |                                  |
|-----------------------|----------------------------------|
| Type                  | aneroid                          |
| Measurement method    | auscultatory                     |
| Measurement range     | 0 - 300 mmHg                     |
| Pressure accuracy     | ± 3 mmHg                         |
| Zero zone             | ± 3 mmHg                         |
| Air release rate      | 2 mmHg/second                    |
| Operating temperature | + 10°C at + 40°C                 |
| Storage temperature   | - 15°C at + 40°C                 |
| Humidity              | 85 % HR                          |
| Weight                | 300 at 400 g according to models |
| Cotton cuffs          | pull on velcro                   |
| Nylon cuffs           | pull on velcro                   |

*Life-time of the Spengler pull on velcro cuff is about 15000 opening/closing operations.*

*In compliance with the standards, the instrument is regulated at the factory according to standard procedures. The end user is not allowed to intervene on the needle position.*

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Hersteller : Spengler S.A.S.